SERIE STEREO

GUIDA UTENTE

INSTRUCTION MANUAL

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
Introduction

The OPTIKA microscopes STEREO series are high performance stereo microscopes with a large possibility of choose between various models with fixed magnification or rotatable objective

Unpacking and Assembling of the Microscope

The components for stereo series are shipped detached for protection. Open the styrofoam protection with care and do not leave any components attached to the packing being removed. Normally this series is shipped ready to be used or you have only to put the eyepieces into their tube. If damage occurs during the transport please contact immediately both the carrier and your supplier.

Alignment and Operation of the Microscope

Interpupillary distance

Move the two eyepiece tubes until only one circular field can be seen through the two eyepieces. If two circles appear the interpupillary distance is too big, if two overlapped circles appear the interpupillary distance is too small.

Focusing the microscope

Try to focus the sample with the magnification 2x by using the focusing knobs or if necessary adjust the high of the microscope along the vertical post.

Magnification

Select the magnification (when is possible) by adjusting the rotatable objective. The total magnification can be calculated by the following equation.

TOTAL MAGNIFICATION = Eyepiece magnification x Objective lens magnification

Stereo head

The binocular tube is built together with the stereo body to form a single piece named "stereo-head".

Stand, illuminator.

The stand is equipped with an inox pole and the focusing mechanism (with stereo head) can be moved vertically to obtain the right focusing.

the fine focusing is obtained by the dedicated knobs.

The STEREO series is formed by models with illumination system and models without it as showed in the summary table. The light can be incident, transmitted and if necessary is possible to use an external cold light source. By using the correct light is possible obtain the better image of your sample.

To take good care of the microscope, the user should avoid the attack of dust and water. If dust and water come to the microscope, sooner or later, fungus will grow and kill the microscope. Please notice that once fungus is grown, even after cleaning, the microscope will have a good chance to be re-infected again. Moreover, grease stain and fingerprints will distort the image, which should be avoided.

Protection against Dust

When the instrument is not used for a long period of time, cover it with the enclosed dust cover. Never leave the eyepiece tube open without any covers. Either leave the eyepiece in the tube (recommended if the microscope is used very frequently) or cover it with wrapping paper or covering cap. Eyepieces and other optical accessories, when not in use, are recommended to be kept in a dry box to protect them against dust and water.

Protection against Water

The instrument should be kept away from water source, pipeline or water sink. Humidity in the room where the instrument located should be as low as possible (relative humidity should be kept below 70%). All the optical accessories are recommended to be kept in a dry box when not in use. The use of dehumidifier and 24-hours air conditioning is highly recommended if the surrounding is too humid.

Cleaning

If dust is found on the optical surface, try to remove by air blower or compress air. For fingerprint, grease stain or dust which cannot be removed by the air blower, 2 possible methods are recommended:

- To breathe lightly on the glass surface and wipe with a clean piece of cloth, lens paper or cotton swab. Please notice that small cotton fibre may be left onto the lens surface if cotton swab is used.
- Use a cotton swab or lens paper, dip with a small amount of absolute alcohol, and clean the lens surface carefully. No other aggressive solvents should be used.

In no circumstances should the user clean any lens surface with dry cotton swab, cloth or lens paper. This will scratch the lens surface causing irreparable damage. Water is not recommended for cleaning of lens as it will leave some water stain on the lens surface and if water residue is left on the lens, fungus can grow causing irreparable damage.

Electrical Parts of the Microscope

Before plugging in the power cord with the supply, make sure that the supplying voltage matches with the operation voltage of the equipment.

Turn off the equipment before plugging in the power cord with the supply.

The best practice is to turn the voltage down to the lowest before turning off the equipment.

Do not turn the power on and off, off and on immediately as this will shorten the life span of the bulb and may cause damage to the electrical system.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.


St-001	Eyepieces Wf5x/22mm
St-002	Eyepieces Wf10x/20 Mm
St-003	Eyepieces Wf15x/15 Mm
St-004	Eyepieces Wf20x/13 Mm
St-005	Micrometric Eyepieces Wf10x/20 Mm
St-040	Dark Field Condenser
St-041	Jewellery Clips
Vc01	Videomicroscopy System
Vc04	Videomicroscopy System With High Resolution

The appliance reports the symbol:

Such symbol means that the appliance can be a precious source of raw materials. Therefore, it must not be disposed of as waste, but separately collected for the recycling and the recovery of the materials contained in it. Such materials, if improperly dispersed into the environment, can be harmful to the environment and to human health.

The producer of the equipment, M.A.D. Apparecchiature Scientifiche – Optika Microscopes, recovers, re-uses and recycles the raw materials contained in the equipment. Such recovery, however, needs your help.

When, at the end of its operating life, you shall decide to dispose of the apparatus, do not try to open it, nor to use parts of it in ways other than reported in this Manual, but bring it back to the Distributor you bought it from, or to other M.A.D. or Optika Microscopes distributors. The Distributor shall collect the apparatus free of charge.

The recovery of the raw materials shall then take place in accordance with the European Directive 2002 / 96 / EC and all other relevant Directives. Never disassemble, nor dispose of as waste, apparatus reporting the “crossed bin” symbol indicated above.

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